Call 0800 102 808 to discuss your specific requirements with one of our product experts.





"It makes getting in and out of bed such a simple matter that I don't really think about it."

Florence, New Plymouth

Phone: 0800 102 808 0800 102 800 Fax: Email: info@multifit.co.nz

Physical address: Multifit International Ltd 10 Pohutu Street Whakatane New Zealand

Postal address: Multifit International Ltd PO Box 748 Whakatane New Zealand

www.multifit.co.nz

MultiFit Multi

Your Independence - Our Innovation

Your Independence - Our Innovation

©Multifit International Ltd 2014

Independence Range

Easy Access Bed A Declaration of Independence



Wherever you place your Easy Access Bed, it's perfect for receiving visitors and enjoying a conversation and a cup of tea or coffee.

Independence Range

Easy Access Bed A Declaration of Independence



The Easy Access Bed sitting, reclining and lying positions

Everyone knows the value of a good nights' sleep. But for people with mobility challenges getting into bed at night and getting up in the morning can be extremely difficult.

Maintaining your mobility is crucial to staying independent and keeping safe at home. If you are recovering from an illness, accident or surgery, or have a chronic health problem, it's important to get mobile and keep as active as possible. This may require assistance from loved ones or care staff, not to mention tremendous effort on the part of yourself. Because many people struggle to get in and out of bed, we thought a great solution would be to create a bed that simply moves you into a semi upright position ready to stand and vice versa.

Multifit is known for innovative products that help people with mobility issues and we are delighted to offer you the Easy Access Bed.

Perfectly tailored to the needs of carers and users alike, it provides a simpler way to get in and out of bed that is easier on you and reduces the need for assistance.

Much more than just a bed or a recliner

Unlike other adjustable beds, the Easy Access Bed is like having many pieces of furniture in one. You save on space and avoid having multiple pieces of furniture and equipment cluttering your home. Control is via a simple hand set with four buttons. Two for up and down control of the whole unit and two more for independent backrest control.

The Easy Access Bed transforms smoothly and seamlessly from a horizontal or a legs elevated bed to a recliner, an upright chair and finally to a chair position that tilts you forwards making it much easier to stand up.




Horizontal bed The easy access bed can be used as a fully horizontal bed, wide enough for you to turn on your side or turn right over for a great night's sleep.

Legs elevated position For those who need to sleep or rest with their legs elevated to improve circulation, a simple manual adjustment allows a 13° elevation of the legrest relative to the seat.



Recliner You can control the angle of recline to obtain the ideal position. Perfect for a nap, watching tv, reading or talking with visitors.

Chair Sit as upright as you wish to have a meal, take medication. or socialize.



Easy to Use

The Easy Access Bed is used just like a normal bed, recliner and chair. Control is via a simple hand set with four buttons, two buttons for up and down and two for independent control of the backrest.



Getting into bed and positioning your bedding is simple

Many people have problems positioning their bedding over their legs and feet once they are in bed. The Easy Access Bed makes this simple. When getting into bed, the starting position is the tilted chair position. Seat yourself as far back as possible so your back is close against the backrest. Notice how the bedding is draped over the handrail. Position your pillow and place the bedding over your lap. Use the controller to recline and the bedding will stay in place as you relax.



Elevating your legs improves circulation The Easy Access Bed offers you the option of reclining or sleeping with your legs elevated.



Independent backrest control

If you like to sit up and keep your legs raised or elevated it is possible to independently control the backrest to obtain the perfect position.



Takes the effort out of getting up

To get up, just do the same in reverse. Move the Easy Access Bed into its tilted chair position, place your feet squarely on the floor, bring your heels back and lean your upper body forward slightly. This puts you in the best position to stand up independently. If you need a little extra help, sturdy armrests can be added which will help you steady yourself and give you something to push up from.



Simple hand set control

The handset is very easy to use with two buttons to control the up and down movement of the whole bed and two buttons to control the backrest independently.

► A a if b





Handy night light

A night light on the side of the bed activated by a switch on the hand set, illuminates the floor if you need to get up at night. Perfect in the bedroom or living room

Who would benefit from the Easy Access Bed?

The Easy Access Bed is an ideal choice if you:

- Have difficulty using a standard bed or a hospital bed.
- Spend day and night in a reclining chair.
- Need to sleep in a 'customised' reclined position or a 'legs elevated' position for health reasons.
- Are recovering from an accident, surgery or illness and you need to get into or out of bed with greater ease.
- Want to stay independent and safe in your own home but this is problematic due to difficulties getting in and out of bed
- Need help transferring from bed to a wheelchair. Because the Easy Access Bed fully supports you into a 'ready to stand' position, it is easier and safer for a carer to assist you from the bed to the chair reducing the risk of injury.





Standard mattress

Deluxe mattress

Pressure relief



Mattress clips

MATTRESSES

The standard mattress is constructed from a base layer of standard foam with a top layer of memory foam for comfort covered with an easy wipe Visco Elastic material.

For a little more luxury choose our deluxe quilted mattress with an additional layer of wool and polyester and an attractive quilted cover. If you require pressure relief talk to us about your specific needs and we can discuss options. Mattresses are easy to remove and replace by unfastening clips that attach the mattress to the frame.



The Easy Access Bed works well in a bedroom or a living room.



ARM RESTS

Optional arm rests provide additional stability and assistance when standing. Wide wooden hand grips provide a stable platform from which you can push yourself up into a standing position. Practical trays are attached to contain items you need close by.



SEAT HEIGHT AND LEG ELEVATION ADJUSTMENT

The Easy Access Bed attaches to its base by means of bolts with four different settings. This enables the seat height to be adjusted to suit your stature and also allows the final (bed) position to be set as either fully horizontal or legs elevated depending on your preference. These adjustments are simple and can be done at home but for your convenience we will pre-set the Easy Access Bed before delivery.

BATTERY BACKUP

A backup battery* can be supplied to protect against power outages by enabling you to get out of bed if there is a power cut.

*A safety back up to get you out of bed if there is a power outage. Not intended to power the Easy Access Bed repeatedly during outages.

Specifications

Easy Access Bed dimensions	1.12m v
Standard frame upholstery	Grey lea
Standard mattress	Memory Visco El
Standard mattress dimensions (W x L x D)	0.9m x 2 50/120F
Deluxe quilted mattress	Memory addition 50/120F
Deluxe quilted mattress dimensions (W x L x D)	0.9m x 2
Pressure relief mattress	Pressur discuss 50/120F
Seat height range	Adjusta
Maximum leg elevation achievable	13° rela
Weight of bed (frame only)	52kg
Weight of mattress	20kg
Maximum user weight	175kg

Easy Access Bed basic unit
Standard mattress with Visco Elastic cover
13° leg elevation adjustment
Independent back rest control
Night light
Deluxe quilted mattress
Pressure relief mattress
Arm rests with side trays (pair)
Battery backup unit

Design changes and improvements may be made without notice.



wide x 2m long when the bed is extended

ather grain Vinyl

ry foam over standard foam with a wipeable lastic cover.

2m x 140mm. HR / 65065V max user weight 125kg

ry foam over standard foam with a quilted cover and nal wool and polyester filling. HR / 65065V, quilted wool top, max user weight 135kg

2m x 160mm.

re relief air mattress options are available call us to your specific needs. HR Roho matress, max user weight 220kg

able from 480 – 600mm

ative to the seat

Standard	Optional
\checkmark	
\checkmark	
\checkmark	
\checkmark	
	\checkmark
	· ·