

Suction rail

| Multifit code | Length |
|---------------|--------|
| MF 3/1X - SK | 300 mm |



Grab Rails

Designed to provide a temporary hand hold the rail simply suctions onto any flat surface. Used to assist with balance in places like the shower bathroom or toilet.

Flip up the release tabs for easy installation and removal.

Easy to attach and relocate.

Note : Not used to support total body weight. May also need to be re-fixed to renew the suction.

Guarantee

- Delivery with in 5 working days or we will shout your whole OT department morning tea.
- The equipment will retain its original function for at least 2 years or we will replace it free of charge.
- Customised & developmental equipment may take longer.

Fixing instructions

- 1/. Flip up the easy release tabs each end of the grab rail.
- 2/. Locate the suction cups over any flat surface.
- 3/. Flip the tabs down and the suction is complete.
- 4/. Check that the rail is secure.

Bed lever - Reinforced

Instruction sheet

Installing the Bed lever



1/. Slide the base of the Bed lever between the mattress and the base of the bed. Position the handle to the most convenient location for your client. This is normally about 600 mm down from the head end of the bed.

Fitting the braided cord & Cleat

2/. Extend the braided cord to the other side of the bed, pull the cord tight through the block and secure through the plastic cleat.

Fitting the strap & buckle


Extend the strap around the base of the bed and join together the clip connector. Pull down on the strap to secure the rail to the bedside.



1
Slim profile against the bed prevents injuries.



2
Cord and cleat tie off on the other side of the bed to prevent the Bed lever from pulling out from the bed.



Parts replacement

| Part Name | Part number |
|-----------------|--------------|
| Steel handrail | MF 3/11 - A |
| Bed lever base | MF 3/11 - B |
| Braided cord | MF 3/1 - BrC |
| Block and cleat | MF 3/1 - BC |

Cleaning

All surfaces are pre finished and only require dusting occasionally. Reissue clean up can be done with any mild detergent and marks can be removed with Jiff.

Demonstrations and training

Available on request. If you would like a demonstration, please phone 0800 102 803

Multifit products are also available for demonstration as required at the Show your ability expo in March each year.

Safety first

Rounded corners prevent injury and snagging on the bed lining.

Slim profile against the bed prevents injuries with no bulky wooden block to catch your legs on when you get out of bed. See figure 1/.

Cord and cleat or strap and buckle fixing prevent the Bed lever pulling out from the bed when going from standing to sitting on the bed.

NOTE : If there is a concern regarding uncontrolled movements in bed then we recommend the "Multifit Safe rail divider" as this complies with resistance to head traps.

Reusability / cost effectiveness

We have tracked the reuse of equipment over the years. Original Multifit raisers made in 1990 are still in good service today.

Contact details

Multifit Hospital Supplies. Ph 0800 102 803—Fax 0800 102 800 - E Mail info@mutifit.co.nz
P O Box 748 - 10 Pohutu, Street, Whakatane

Mission statement

Quality of life is precious. Our mission is :

To provide you with useful, quality, and long life equipment, that enables your clients to make the most of their lives.

Serving people with disabilities since 1990