Double bed rail for wire wove bed MF3/3

Multifit Hospital Supplies . Enquiries Ph 07 3080045 8 am – 3 pm, Fax 0800 102 800, E Mail multifit@wave.co.nz

Specially designed to clamp & tie to the sides of wire wove type beds, or sectional hospital beds.

Adjustable in width

These rails are connected with a telescopic section that enables you to actually clamp it to the sides of the bed frame, making it easy for you to get a customized fit for your client every time.

Mattress comfort remains.

Because the clamping rail extends under the mattress, it can not be felt through the mattress by your client

Stability up & down , side to side is solid

By clamping to the bed side rails, the rail takes on the rigidity of the bed itself. The clamping rail that extends under the bed is designed to be self bracing, therefor providing rigidity of the rail side to side

Requires no tools to fit

You don't need to drill any holes or fit any bolts, in fact you done even need any tools to fit it. Just turn the knob to adjust the width and activate the clamping action, then tie down the straps and its done.

Reusable because its adjustable and because your clients bed is protected by rubber liners where the rail touches the bed side your clients bed will not be damaged by having the rail fitted





Here the rail is fitted to a wire wove bed that has 35 mm x 35 mm steel angle iron side frames. Being adjustable in width the rail will fit most width single beds

Some wire wove beds have a timber side frame. The clamping of the rail works just the same, but you extend the tie straps By turning this knob the rail is adjusted in



Telescopic clamping rail extends under the bed to provide rigidity for the rail and comfort for your client

Specifications :

- Adjusts from 1150 mm to 760 mm wide so you can get that customized fit . Made to measure rails can also be made.
- The height of the rail extends 470 mm above the bed frame, and is 225 mm wide, & 25 mm diameter, so that your client gets a plenty of rail to grasp while not being too bulky at the bed side..
- Made from 1.6 mm steel for strength, mig welded for strong joint strength , etch primed & powder coated for long lasting good looks and hygiene.
- Fitting of the rail is by the adjustable knob and the braided strap, which means there are no tools required to fit the rail.
- Where the rail attaches to the bed we provide a rubber liner so that your clients bed is not damaged.



The rail fitted to a clients bed

Made to order only