

Bed Block MF 0/1 & MF 0/1C



Multifit

Your Independence - Our Innovation

Fitting Instructions

PLEASE READ - before installation. If in doubt over the fixing of this Raiser, please contact Multifit directly.

CHOOSE THE CORRECT RAISER - As chair, sofa and bed designs are many and varied, it is important to choose a Raiser that best suits the type you have. Please refer to our website to see chair, sofa and bed types and the Raisers that fit them. We are also able to custom make raisers to suit, and fit then on site if required.



USER PRECAUTIONS - This product is only to be installed for the purpose for which it is intended (for furniture placed on a flat, stable and horizontal surface). If your raised furniture sits on a polished floor and you are concerned about it slipping around on the floor, then please order and fit our "non slip adhesive foot pads" to the bottom of the raiser. This raiser is for indoor use only.

MAINTENANCE - Check tightness of wing nuts / knobs and fixings periodically.

CLEAN - with damp cloth. Use non-abrasive detergent and dry thoroughly.

MOVING - Only move the raised furniture when it is unoccupied. After moving check that the raised furniture remains fitted correctly.

ASSESSMENTS - Clients who are very physical in the way they use their furniture can impose very high loads on the raised furniture. Part of the risk assessment should take this interaction into account when choosing the appropriate raiser.

MAXIMUM USER WEIGHT - 380 kg including the furniture. Please contact us for higher user weight limits.

DEMONSTRATIONS - available on request phone 0800 102 803 (includes the Show your ability expo).




PARTS REPLACEMENT - Refer to www.multifit.co.nz for a complete parts list.

OPTIONAL EXTRAS - If the legs of the bed are larger than the opening in the Bed Block then use our Wide pad Bed Raiser. If the bed is a divan type bed then use our Divan type Bed Raiser. See www.multifit.co.nz


MISSION - Our mission is to provide you with the safest, most cost effective raiser to meet your needs.

Fitting Bed Blocks Safely

- 1 Lift the foot end of the bed with the Multifit **Bedjack** and place a Bed block under the feet in the middle of the bed.



- 2 Remove the Multifit **Bedjack**. Place a weight at the foot end of the bed. This will cause the head end of the bed to rise. Place bed blocks under the legs at the head of the bed.



- 3 Lower the head end & place a raiser at the foot. Then using the **Bedjack** raise the head end and place the 3rd raiser

