

# Single & Double Bed Raisers



## MF 2/4 & 2/2 Fitting Instructions

**Multifit**  
Your Independence - Our Innovation

**PLEASE READ** - before installation. If in doubt over the fixing of this Raiser, please contact Multifit directly.

**CHOOSE THE CORRECT RAISER** - As chair, sofa and bed designs are many and varied, it is important to choose a Raiser that best suits the type you have. Please refer to our website to see chair, sofa and bed types and the Raisers that fit them. We are also able to custom make raisers to suit, and fit then on site if required.



**USER PRECAUTIONS** - This product is only to be installed for the purpose for which it is intended (for furniture placed on a flat, stable and horizontal surface). If your raised furniture sits on a polished floor and you are concerned about it slipping around on the floor, then please order and fit our "non slip adhesive foot pads" to the bottom of the raiser. This raiser is for indoor use only.

**MAINTENANCE** - Check tightness of wing nuts / knobs and fixings periodically.

**CLEAN** - with damp cloth. Use non-abrasive detergent and dry thoroughly

**MOVING** - Only move the raised furniture when it is unoccupied. After moving check that the raised furniture remains fitted correctly.

**ASSESSMENTS** - Clients who are very physical in the way they use their furniture can impose very high loads on the raised furniture. Part of the risk assessment should take this interaction into account when choosing the appropriate raiser.

**MAXIMUM USER WEIGHT** - 380 kg including the furniture. Please contact us for higher user weight limits

**DEMONSTRATIONS** - available on request phone 0800 102 803 (includes the Show your ability expo)



**PARTS REPLACEMENT** - Refer to [www.multifit.co.nz](http://www.multifit.co.nz) for a complete parts list.


**OPTIONAL EXTRAS** - If the legs of the bed are larger than the opening in the Raiser then use our Wide pad Bed Raiser. If the bed is a divan type bed then use our Divan type Bed Raiser. See [www.multifit.co.nz](http://www.multifit.co.nz)

**MISSION** - Our mission is to provide you with the safest most cost effective raiser to meet your needs.


## Fitting Bed Raisers Safely

**1** At the foot of the bed, adjust the raiser out so the pads of the raiser align with the center of each foot of the bed. Tighten the wing nuts. This part of the process is the same for both Single and Double bed raisers.

**2** Beds with feet at the foot & head end use this method.




**3** Raise the head end of the bed from the side of the bed. Place a weight on the near foot end corner. Place a raiser under the feet at the head end. Lower the Jack.




**4** Raiser the foot end of the bed. Fit the raiser. Lower the Jack. Check that the bed feels stable and safe. **This raiser will also fit bed ends as pictured to the right.**




**2** Beds with a center row of feet use this method.



**3** With the Multifit Bed Jack raise the end of the bed. Place the Bed Raiser under the center legs. Lower and remove the Jack.



**4** Put a weight on the foot end of the bed (anything 3-5kg will do). Place a raiser under the head end. Remove the weight.



**5** Raise the end of the bed again. Fit the last raiser. Lower the bed. Check that the bed feels stable and safe.

