



## Fitting Instructions

**PLEASE READ** - before installation. If in doubt over the fixing of this Raiser, please contact Multifit directly.

**CHOOSE THE CORRECT RAISER** - As chair, sofa and bed designs are many and varied, it is important to choose a Raiser that best suits the type you have. Please refer to our website to see chair, sofa and bed types and the Raisers that fit them. We are also able to custom make raisers to suit, and fit then on site if required.



**USER PRECAUTIONS** - This product is only to be installed for the purpose for which it is intended (for furniture placed on a flat, stable and horizontal surface). If your raised furniture sits on a polished floor and you are concerned about it slipping around on the floor, then please order and fit our "non slip adhesive foot pads " to the bottom of the raiser. This raiser is for indoor use only.

**MAINTENANCE** - Check tightness of wing nuts / knobs and fixings periodically.

**CLEAN** - with damp cloth. Use non-abrasive detergent and dry thoroughly.

**MOVING** - Only move the raised furniture when it is unoccupied. After moving check that the raised furniture remains fitted correctly.

**ASSESSMENTS** - Clients who are very physical in the way they use their furniture can impose very high loads on the raised furniture. Part of the risk assessment should take this interaction into account when choosing the appropriate raiser.

**MAXIMUM USER WEIGHT** - 380 kg including the furniture. Please contact us for higher user weight limits.

**DEMONSTRATIONS** - available on request. Phone 0800 102 803 (includes the Show your ability expo).



**PARTS REPLACEMENT** - Refer to [www.multifit.co.nz](http://www.multifit.co.nz) for a complete parts list.

**OPTIONAL EXTRAS** - If the legs of the Sofa are wider than the opening in the Raiser then use our Wide pad Sofa Raiser. If the Sofa is a divan type then use our Divan type Sofa Raiser. See [www.multifit.co.nz](http://www.multifit.co.nz)

**MISSION** - Our mission is to provide you with the safest, most cost effective raiser to meet your needs.

## Fitting Sofa Raisers Safely



**WARNING Lifting Hazard** - Single person lift could cause back injury. We recommend using the Multifit Bed Jack. This Jack is designed to work with the Multifit line of raisers in that the space under the telescopic section of each raiser allows the foot of the jack to be placed.

**2** With the Multifit Bed Jack raise the end of the Sofa. Place the Sofa Raiser under the feet of the Sofa. Lower and remove the Jack.

**3** Raise the other end of the Sofa. Fit the last raiser. Lower the Sofa. Check that the Sofa feels stable and safe

**1**

At the end of the Sofa, adjust the raiser out so the pads of the raiser align with the each foot. Tighten the wing nuts.

