

Call 0800 102 803 to discuss  
your specific requirements with  
one of our product experts.

**The Multifit Guarantee covers all our Chair and Bed  
Raisers**

Delivery within 5 working days from order or we will  
treat your entire department to a morning tea.

The equipment will fit or we will customise one  
that will.

The equipment will retain its original function for at  
least two years or we will replace it at no extra cost.

Phone: **0800 102 803**

Email: **reception@multifit.co.nz**

Postal address:

**Multifit Hospital Supplies Ltd,**

**PO Box 748**


**Whakatane**

**New Zealand**

**[www.multifit.co.nz](http://www.multifit.co.nz)**



©Multifit Hospital Supplies Ltd 2019



**Multifit**

Your Independence - Our Innovation

Phone 0800 102 803

[www.multifit.co.nz](http://www.multifit.co.nz) [reception@multifit.co.nz](mailto:reception@multifit.co.nz)




# Multifit Catalogue

Version: 1



Innovative solutions that  
help you get in and out of  
a bed or chair



## Multifit

Your Independence - Our Innovation

Phone 0800 102 803

[www.multifit.co.nz](http://www.multifit.co.nz) [reception@multifit.co.nz](mailto:reception@multifit.co.nz)

# Contents

---

## **RAISERS**

LazyBoy Chair Raiser .....	3
Standard Chair Raiser .....	4
Big Softie Chair Raiser .....	5
Sofa Raisers .....	6
Widened Raisers.....	6
Customized Raisers.....	6
Bed Raisers - Single and Double.....	7
Bed Blocks.....	8
Divan Raiser.....	8
X Chair Raiser.....	9
Five Leg Chair Raiser.....	9
Pedestal Chair Raiser .....	10
Open Front Chair Raiser .....	10

## **LEVERS AND RAILS**

Bed Lever.....	11
Chair Side Transfer Rail .....	11
Ensuring Safety with a Bed Lever.....	12
Bed Side Transfer Rail .....	13
Profiling Bed Transfer Rail.....	13

## **MISCELLANEOUS**

Assessment Handrail.....	14
Handybar .....	14
Shower Remote .....	14
Slide Sheet .....	15
Bed Jack.....	15

## **MULTIFIT INTERNATIONAL**

LegLifter Transfer Rail.....	16
LegLifter .....	17

# About Multifit

---

First established in 1992, Multifit is a New Zealand owned and run business that has pioneered the design and development of many new and life changing assistive solutions and devices. Products are patented both in New Zealand and internationally and Multifit has been exporting since 2005. Director and Head of Product Development, Wayne Manson, is well attuned to the needs of end users and clinical staff including Assessors, Occupational Therapists and Rest Home and Retirement Village Clinical Managers.

*“Multifit has been enhancing independence and quality of life for the elderly and disabled through our innovative mobility solutions for over 25 years and this is something we are very proud of. Our Chair Raisers are safe, robust, durable and attractive. They are cost effective products that stand the test of time and are frequently used many times over. We also understand that choosing the appropriate Chair Raiser ensures the most cost effective solution to raising a chair in the safest possible way.”*



**Wayne Manson**

*Director*



# Multifit Hospital Supplies

## *Innovative solutions that help you get in and out of a chair or bed*

The Multifit Hospital Supplies range of sofa, chair and bed raisers, levers and rails has been developed to assist mobility and promote independence for the end user as well as supporting the care professional. Fully adjustable to ensure a perfect fit, Multifit Chair Raisers are made in New Zealand from furniture grade timber and colour coated sheet metal for long lasting good looks. All our raisers adjust by hand with no tools required and prefinished surfaces make cleaning easy and hygienic.

Safety and security are paramount in all our designs and attention to detail has gone into each product to ensure maximum stability in terms of attaching the product securely to furniture to eliminate undesired movement.

A problem that many Therapists will recognise is arriving at a client's home only to find that the product they have brought does not fit the chair or bed. Our signature "telescopic" design overcomes this common problem making our chair and bed raisers fully adjustable over a wide range of furniture sizes.

The robust and solid manufacture of our products make it the economical choice. Our products are so durable they frequently outlast their owners and can be re-used over and over.

### **Simple tips to make getting up from your chair easier**

Many of us have difficulty standing up from chairs as we get older. Here are some simple tips to that may help you.

- Older people have problems standing from a seat less than 475mm from the floor, so try sitting in a chair with a higher seat or a firmer seating cushion.
- Solid arm rests that extend forward of the seat are best. Ideally you should be able to grip around the end of the arm rest with your fingers. This will help with stability as you try to stand and give you something to push up from with your arms.
- Good technique is important. Move as far forward to the edge of the seat as you can and tuck your heels back as close to the front of the chair as possible. Spacing your feet apart will help with stability.
- Tilt your upper body forward a little – a good idea is to say to yourself 'nose over toes' as you push up from your heels to stand.

- A person's size relative to the chair is also important for proper support. For example, a smaller person would struggle to stand up from a wide chair with widely spaced arm rests. In this instance select a narrower chair. For a heavier person, a chair with a deeper and wider seat and a more upright back may be better.

**If you have explored various options with your existing furniture then raising your seat height using a Multifit Raiser is a very cost effective way to improve your ability to stand independently.**

### **Rails and handholds can make a big difference**

If your seat is the correct height and you are using a good technique but are still struggling to stand, think about how you are using your arm rests. The best arm rest extends beyond the front of the chair, is comfortable to hold and is firm to push up from.



There may be other chairs in the room which have better arm rests to push up from. Sometimes minor changes in arm rest location can make a significant improvement to your ability to push up from, or maintain balance when getting on and off a chair.

Arm rests which are widely spaced are difficult to push up from, so look for chairs which are narrower and have arm rests which are close to the sides of your body as you stand.

Multifit Chair Transfer Rails are a handy solution that provide that extra support and security. These simply slide under the chair base and the hand holds can be adjusted to suit your requirements.

#### Calculating the correct height for a chair raiser

- Calculate the correct seat height by measuring the distance from the floor to the crease at the back of your knees. Then compare this measurement with the seat height of your chair.
- When seated, your hips and knees should be at right angles with your feet flat on the floor wearing normal indoor shoes or slippers.
- A slightly higher seat can make it easier to stand up and sit down, particularly if you have any pain or weakness in your legs.



Correct height



Too low



Too high

# Chair Raisers



KEY BENEFITS

Safe and secure  
Long lasting  
Robust, re-usable, economical  
Telescopic adjustability for perfect fit  
Unique Multifit Tri-knob will not loosen

## LAZY BOY CHAIR RAISER

This chair raiser suits Lazy Boy type rocking chairs which have a wooden base frame similar to the one pictured. The raiser is fully adjustable to fit a variety of base sizes. Because the raiser has a special device that locks it onto the chair base, the chair has perfect stability but the user can still rock and use the reclining facility.



### Special features for added safety

Safety and stability are especially important with rocking chairs as the rocking and reclining action can cause the chair to come off a standard raiser.

- Splayed feet widen the raiser base to increase stability
- Strong clamps at the rear secure the raiser tightly to the chair. When the chair tilts forward, the clamps ensure that the base of the chair stays on your raiser and eliminates the risk of tipping.
- To enhance safety and security this product utilises our unique Multifit Tri-knob, manufactured exclusively for our products and designed with feedback from New Zealand Occupational Therapists. Due to its unique shape, the Tri-knob enables maximum grip and leverage ensuring it does not loosen over time. It has rounded edges and a slim profile to ensure that users will not snag themselves as they walk past.



### Installation

1. Loosen the four wing nuts on the sides of the raiser so you can extend it to a larger size than needed. Tip the chair forwards so the base is visible and adjust the raiser to fit over the feet of the chair. Tighten the wing nuts firmly.
2. With the spreader bar placed over the timber rockers of the chair, tighten the clamp screws to secure the spreader bar down onto the rockers by turning the top wing nut and then the lock nut below.
3. Tip the chair back onto its feet and rock it to make sure it is safe and secure before using.



### Specifications

Min max base adjustment	490 x 528mm / 705 x 880mm
Heights available	50, 75, 100 and 150mm
Maximum user weight	380kg
Standards	Based on BS EN 1022: 1997 stability of furniture.

### Order Codes

Code	Height
MF 1/3 - 50	50mm
MF 1/3 - 75	75mm
MF 1/3 - 100	100mm
MF 1/3 - 150	150mm



# Chair Raisers



## KEY BENEFITS

Safe and secure  
Long lasting  
Robust, re-usable, economical  
Telescopic adjustability for perfect fit  
Unique Multifit Tri-knob will not loosen

## STANDARD CHAIR RAISER



For square shaped lounge chairs that do not rock or recline, the Standard Chair Raiser is a perfect solution. Providing a robust, safe platform to raise the chair, the raiser adjusts smoothly to fit the dimensions of the chair base and wing nuts clamp tightly to fix the raiser at the size required.

### Special features for added safety

To enhance safety and security this product utilises our unique Multifit Tri-knob, manufactured exclusively for our products and designed with feedback from New Zealand Occupational Therapists. Due to its unique shape, the Tri-knob enables maximum grip and leverage ensuring it does not loosen over time. It has rounded edges and a slim profile to ensure that users will not snag themselves as they walk past.

### Installation

1. Loosen the four wing nuts on the sides of the raiser.
2. Tip the chair forward so the base of the chair is upturned and adjust the raiser to fit over the base of the chair.
3. Tighten the wing nuts on the raiser firmly.
4. Place the raiser on the floor behind the chair and slide the chair back onto the raiser so that the base of the chair fits within the rim of the raiser. Make sure the chair feels stable before using.



### Specifications

Min max base adjustment	Standard: 528 x 483mm / 873 x 816mm (MF 1/1) Mini: 385 x 435mm / 500 x 645mm (MF 1/0) Large : 680 x 628mm / 1020 x 956mm (MF 1/2)
Heights available	50, 75, 100 and 150mm
Maximum user weight	380kg
Standards	Based on BS EN 1022: 1997 stability of furniture.

### Order Codes

Code	Height
MF 1/1 - 50	50mm
MF 1/1 - 75	75mm
MF 1/1 - 100	100mm
MF 1/1 - 150	150mm

# Chair Raisers



## KEY BENEFITS

Safe and secure  
Long lasting  
Robust, re-usable, economical  
Telescopic adjustability for perfect fit  
Unique Multifit Tri-knob will not loosen

## BIG SOFTIE CHAIR RAISER

This raiser suits Big Softie type chairs which have a tubular steel base frame. The raiser adjusts smoothly to fit the dimensions of the base of your chair and fixes securely. Providing you with a robust, safe platform to raise the chair, our special cord and cleat fixing device ensures that the raiser is safely attached and that the user can still rock and recline with confidence







### Special features for added safety

Safety and stability is especially important with rocking chairs as the rocking and reclining action can cause the chair to come off a standard raiser.

- Splayed feet widen the raiser base to increase stability
- A secure cord and cleat fixing at the rear of the raiser secures it tightly to the chair base frame. When the chair tilts forward, this fastening ensures that the base of the chair stays on your raiser and eliminates the risk of tipping.
- To enhance safety and security this product utilises our unique Multifit Tri-knob, manufactured exclusively for our products and designed with feedback from New Zealand Occupational Therapists. Due to its unique shape, the Tri-knob enables maximum grip and leverage ensuring it does not loosen over time. It has rounded edges and a slim profile to ensure that users will not snag themselves as they walk past.

### Installation

1. Loosen the four wing nuts on the sides of the raiser.
2. Roll the chair forward so the base is upturned and adjust the raiser to fit over the feet of the chair. Tighten the wing nuts firmly.
3. Tie the raiser to the chair base using the cord and cleats as follows:
  - Go over then under the steel tube A
  - Go over then around the chair raiser rail B
  - Thread the cord through the hole in cleat C
  - Finish by pulling the cord upwards into the wedge D
4. Roll the chair back onto its feet and test that it is safe and stable before using.



### Specifications

Min / max base adjustment	490 x 528mm / 705 x 880mm
Heights available	50, 75, 100 and 150mm
Maximum user weight	380kg
Standards	Based on BS EN 1022: 1997 stability of furniture.

### Order Codes

Code	Height
MF 1/7 - 50	50mm
MF 1/7 - 75	75mm
MF 1/7 - 100	100mm
MF 1/7 - 150	150mm

# Chair Raisers



## KEY BENEFITS




Safe and secure  
Long lasting  
Robust, re-usable, economical  
Telescopic adjustability for perfect fit  
Unique Multifit Tri-knob will not loosen

## SOFA RAISER

Sofa Raisers are blocks that are connected with a telescopic midsection to enable adjustment to the particular size of the sofa. The base area of each block has a 20mm high rim and rubber surface to eliminate any movement of the sofa foot. Splayed feet on the blocks increase safety and stability.

### Special features for added safety

To enhance safety and security this product utilises our unique Multifit Tri-knob, manufactured exclusively for our products and designed with feedback from New Zealand Occupational Therapists. Due to its unique shape, the Tri-knob enables maximum grip and leverage ensuring it does not loosen over time. It has rounded edges and a slim profile to ensure that users will not snag themselves as they walk past.



### Specifications

Adjustment range	624 x 900mm
Heights available	50, 75, 100 and 150mm
Maximum user weight	380kg
Non slip rubber pads	90 x 90mm
Lip around the foot pad	20mm high
Standards	Based on BS EN 1022: 1997 stability of furniture.

### Order Codes


Code	Height
MF 2/1 - 50	50mm
MF 2/1 - 75	75mm
MF 2/1 - 100	100mm
MF 2/1 - 150	150mm

### Widened Raisers

Multifit Sofa, Bed and Divan Raisers can be adjusted per request to fit extra wide feet. Raisers can only be widened for heights 75 mm and above.

### Customized Raisers

If you have a chair, bed or sofa .. or anything else, that you would like to raise, we can make a raiser to suit your needs! Quotes on request.





# Bed Raisers



KEY BENEFITS

Safe and secure  
Long lasting  
Robust, re-usable, economical  
Telescopic adjustability for perfect fit  
Unique Multifit Tri-knob will not loosen

## BED RAISERS SINGLE & DOUBLE




Multifit Bed Raisers come in single and double sizing and are easier and quicker to fit than traditional bed blocks with the added benefit of greater safety and stability.

### Special features for added safety


To enhance safety and security this product utilises our unique Multifit Tri-knob, manufactured exclusively for our products and designed with feedback from New Zealand Occupational Therapists. Due to its unique shape, the Tri-knob enables maximum grip and leverage ensuring it does not loosen over time. It has rounded edges and a slim profile to ensure that users will not snag themselves as they walk past.

### Installation

1. Adjust the raisers to suit the overall width of the bed feet and tighten the wing nuts.
2. Lift the foot end of the bed. If you are fitting the raiser to a double bed with a centre set of feet its easier to fit the raiser to the centre set of feet first.
3. With the centre raiser in place, lower the head end of the bed and place a raiser under the feet at the head end of the bed. Finally lift the foot end and place the final raiser.



**NOTE:** To eliminate the risk of manual handling injury when lifting a bed to fit raisers, please use our Multifit Bed Jack. It's simple to use, will easily lift a bed in minutes and will ensure your safety (see page 14)



### Specifications

Adjustment range	Single Bed: 824 x 1100mm Double Bed: 1150 x 1700mm
Heights available	50, 75, 100 and 150mm
Maximum user weight	380kg
Non slip rubber pads	90 x 90mm
Lip around the foot pad	20mm high
Standards	Based on BS EN 1022: 1997 stability of furniture.

### Order Codes

Type	Height	Set of 2	Set of 3
Single Double	50mm	MF 2/2-50/2 MF 2/4-50/2	MF 2/2-50/3 MF 2/4-50/3
Single Double	75mm	MF 2/2-75/2 MF 2/4-75/2	MF 2/2-75/3 MF 2/4-75/3
Single Double	100mm	MF 2/2-100/2 MF 2/4-100/2	MF 2/2-100/3 MF 2/4-100/3
Single Double	150mm	MF 2/2-150/2 MF 2/4-150/2	MF 2/2-150/3 MF 2/4-150/3

# Bed Raisers



**KEY BENEFITS**  
Safe and secure  
Long lasting  
Robust, re-usable, economical

## BED BLOCKS

Lightweight Multifit Bed Blocks are pyramid shaped to provide maximum stability with a convenient finger hole to facilitate carrying.



### Specifications

Recess for bed feet	75 x 75 x 40mm deep
Maximum user weight	380kg over 4 feet

### Order Codes

Code	Height
MF 0/1	50mm
MF 0/2	75mm
MF 0/3	100mm
MF 0/4	150mm

# Special Chair Raisers

**KEY BENEFITS**  
Safe and secure  
Long lasting  
Robust, re-usable, economical  
Telescopic adjustability for perfect fit  
Unique Multifit Tri-knob will not loosen

## DIVAN RAISER

The Divan Raiser is designed for beds or sofas that have no feet, gliders or castors. Each foot pad has two rims that allow the sofa or bed corner to sit snugly onto the foot pad. Self-bracing blocks are linked with a telescopic section that connects with a cord.



### Specifications

Divan raiser type	Code	Heights
Sofa	MF 2/5	50, 75, 100 or 150mm
Single Bed	MF 2/2 D	50, 75, 100 or 150mm
Double Bed	MF 2/4 D	50, 75, 100 or 150mm



# Special Chair Raisers

### KEY BENEFITS

Safe and secure  
Long lasting  
Robust, re-usable, economical  
Telescopic adjustability for perfect fit  
Unique Multifit Tri-knob will not loosen

## X CHAIR RAISER

Fully adjustable to ensure a perfect fit. Made from furniture grade timber and colour coated sheet metal for long lasting good looks. Adjusts by hand, no tools required. Prefinished surfaces make cleaning easy and hygienic.



### Specifications

Heights	50, 75, 100 or 150mm
Non-slip rubber pad	90 x 90mm
Radius adjustment range	330 x 420mm
Size adjustment range	645 x 645mm - 900 x 900mm
Order Code	MF 1/X




## 5 LEG CHAIR RAISER

Fits a chair with a five pronged base. Fully adjustable for a perfect fit.



### Specifications

Heights	50, 75, 100 or 150mm
Non slip rubber pad	90mm x 90mm
Radius adjustment range	330mm x 420mm
Order Code	MF 1/Q





# Special Chair Raisers

**KEY BENEFITS**  
Safe and secure  
Long lasting  
Robust, re-usable, economical  
Telescopic adjustability for perfect fit  
Unique Multifit Tri-knob will not loosen

## PEDESTAL CHAIR RAISER

Designed to fit chairs with a round base. Diameter adjusts so you get a perfect fit every time. Splayed feet increase stability. Made from furniture grade timber and colour coated sheet metal for long lasting good looks and easy cleaning. The brackets can be modified to fit Stressless Recliners from Danskemobler.



### Specifications

Heights	50, 75, 100 or 150mm
Pad clips	20mm high
Radius adjustment range	530 x 620mm (700mm if you remove then reinstall the rim clips)
Order Code	MF 1/PED



These are custom made for each individual pedestal chair, due to the number of variations in chairs available. Before supplying, we need the following from you. 1) Photo of the base - at ground level. Place chair on a hard surface so we can see what's underneath the chair. 2) The Diameter of the base 3) Height of the edge of the base (from the ground) - where the clips are to fit. We need to put the clips on the fixed part of the base so we don't stop the chair swivelling, and this varies depending on the type of pedestal chair.

## OPEN FRONT CHAIR RAISER

The front of this raiser is recessed enabling you to get your feet further under the front edge of the chair. This means that you do not have to slide as far forward on the seat in order to get your body weight over your feet to stand and you can also make better use of the arm rests.



### Specifications

Heights	50, 75, 100 or 150mm
Lip at the corners	20mm high
Width adjustment range	500 x 700mm
Length adjustment range	620 x 905mm
Order Code	MF 1/OF



# Bed Levers & Rails

**KEY BENEFITS**  
 Safe and secure  
 Long lasting  
 Robust, re-usable, economical

## BED LEVER

Placed between the mattress and the base of your bed, the Bed Lever provides a useful handhold to assist you on and off your bed safely. Correct fitting is essential to the safe use of this product. Take care to read the fitting instructions and ensure that it is refitted the same way every time the bed is made.



**Bed Lever Wedge**  
*Has 30% better lateral stability*

Base size	460x800x16mm
Water resistant compressed sheet	
Handle height	560mm
Max user weight	150kg
Order code	MF 3/1AS-CC



**Bed Lever Oblong**

Base size	340x800x16mm
Water resistant compressed sheet	
Handle height	560mm
Max user weight	100kg
Order code	MF 3/11S-CC



Fig 1 - Position the Handle at about hip level when lying. Do not place the handle near the pillow.

Fig 2 - Fit the Handle tightly against the side of the mattress by firmly connecting the Cord from the Block, across and under the mattress to the Plastic Cleat on the Bed Lever base. For increased Handle stability add the Stabilizing Support (opposite page).



## CHAIR SIDE TRANSFER RAIL

Chair Side Transfer Rails assist with independent access in and out of your chair. The height of the rail handles can be adjusted to suit your stature and the base of the rail extends under the chair for stability. Right and left hand and double sided versions are available.

**Order Codes**

Item	Code
Hand rail both sides	MF 1/9 D
Left hand side	MF 1/8 L
Right hand side	MF 1/8 R



# Ensuring Safety With a Bed Lever

**Preventing head traps** - If you've made your assessment and think a Bed Lever is appropriate, then your number one concern should be to stop the handle pulling out from the bed where the person might roll out of bed, pushing the handle out enough to become trapped in the gap. *See Fig 1.* So make sure it's fixed securely and you can not force your forearm between the handle and the side of the mattress. When positioning the handle keep it at least 318 mm down from the end of the bed and well away from the pillow.

**Does the handle wobble too much?** - The standard Bed Lever with the Cord and Cleat suits people who are using the handle to assist with balance. *See Fig 2.* But if you think the Bed Lever handle wobbles too much or if you are at all concerned that your client might push or pull on the handle with a degree of force, then adding the "Bed Lever Stabilising support, MF - CCBFX" will anchor the Bed Lever in place, greatly improving the handle stability and more safely secure the handle tightly to the bedside. *See Fig 3.* This is a strong plastic moulding that straps to the base of the bed and bolts to the base of the Bed Lever. The Cord and Cleat are replaced with a tubular steel mattress Back Fixing which holds tightly to the mattress. This method offers the highest degree of Bed Lever safety while still allows the Bed Lever to be loosened so the bed can be made.

**Refitting** - If at any time the Bed Lever is removed or loosened from the bedside, take care that it is re fitted to the original standard before using.

**Consider the user's condition** - 1) It is common for people to move in their sleep. Some move very little, while others move a lot or have uncontrolled movements. 2) Often people rely on the Bed Lever handle to transfer on and off their bed in the dark of night. Take these things into consideration when choosing a Bed Lever. While Bed Levers can be very helpful, they might not be suitable for everyone. People who have a high degree of violent, or uncontrolled movements, mental or physical issues, might be better suited to sleeping in an Adjustable Medical bed.

**Retrofitting old bed lever handles** - Adding an infill to our older style Bed lever handles will resist head traps. *See Fig 4.*




Fig 1




Fig 2 - Does the handle wobble too much?




Fig 3  
Stabilising Support  
MF-CCBFX

Bed Lever base bolts to the plastic moulding

Straps hold the plastic moulding to the base of the bed

Fig 4  
Bed Lever handle infill .  
Code MF 3/11 - IF




Fig 4



# Bed Levers & Rails

**KEY BENEFITS**  
Safe and secure  
Long lasting  
Robust, re-usable, economical

## BED SIDE TRANSFER RAIL

Bed Side Transfer Rails assist with rolling over in bed, sitting up and getting in and out of bed. A handrail extends outwards at right angles to the bed which is useful for transferring from a wheelchair or walking frame. Options for right or left hand sides of the bed are available.

### Specifications

Rail size	810h x 600w x 32mm diameter
Max user weight	130kg

### Order Codes

Right hand	MF BTR/R
Left hand	MF BTR/L



## PROFILING BED TRANSFER RAIL

The Profiling Bed Transfer Rail is made to fit a profiling bed including the popular Etude Scan bed. The rail extends perpendicular from the bedside to provide a handhold in front of the seated person. Fixes with Allen key bolts.

### Specifications

Extends out from bedside	300mm
Rail diameter	25mm
Rail height	530mm
Max user weight	130kg
Order code	MF ETR



# Miscellaneous

## ASSESSMENT HANDRAIL

Not intended for the support of a person, but ideal to assess where a permanent rail should be fixed as different angles and positions can be tried.

Order code	MF 11/1
------------	---------



## HANDYBAR

Create a safe and secure handhold to assist you when getting in or out of your car. The Handy Bar simply inserts into the door latch U-shaped striker on the car door frame.

Order code	MF HB
------------	-------



## SHOWER REMOTE

An adaptation of the Felton type shower mixer lever that enables the shower to be turned on and off from a lower position. Ideal for those with arm or shoulder problems. The Shower Remote will fit a standard shower cubicle only, as the remote needs to be screwed into the wall or timber framing. The advantage is that the water can be turned on and off using one hand only with an accurate adjustment of the hot /cold water.



Installation	Builder or plumber installs
Easy of use	Can be use with one hand
Order code	MF 6/3A

Installation	Anyone can install
Easy of use	Requires two hands to use
Order code	MF 6/3B

# Miscellaneous

## SLIDE SHEET

Slide sheets are used to facilitate the movement or transfer of patients on beds, chairs and in showers. They assist staff to practice 'No Lift' manual handling procedures in line with Work Cover, Occupational Health and Safety policies.

### Order Codes

1m x 1.5m	MF SS - AMS 1
2m x 1.5m	MF SS - AMS 2



## BED JACK

The Multifit Bed Jack is ideal for reducing manual handling risks when fitting bed raisers. Simply place the foot of the Jack under the edge of the bed and push the button to effortlessly raise or lower the bed.

Lift stroke	20 x 600mm
Width	210 x 870mm
Lifting capacity	250kg
Weight of jack	7kg





# Multifit International

## LEG LIFTER TRANSFER RAIL

A bedside transfer rail to assist you sitting and standing up from your bed is available made specifically for use with the Multifit LegLifter.

### Specifications and Order Codes

Extends out from the bed	276mm
Rail height from mattress	427mm
<b>Order codes</b>	
Right side of bed	MK3LLTRL
Left side of bed	MK3LLTRR




Leg Lifter fitted to an innerspring bed.



## MK3L LEG LIFTER FITTED TO A PROFILE BED

With the addition of profile bed clamps, the Leg Lifter can be fitted to the side rail of most profile beds. The standard clamp will fit bed side rails up to 30mm wide x 60mm deep.



Profile bed clamp



Leg Lifter fitted to the side rail of a profile bed.

# LEG LIFTER

The Leg Lifter is a power-driven appliance designed to gently lift a person’s legs from a bedside seated position to the surface of the bed. It has proven to be a successful mobility solution in providing independence for elderly and disabled clients who have specific difficulty getting their legs into bed. For therapists and caregivers, it saves considerable time and physical exertion.

It can be fitted to some profiling beds and most domestic bed types, including innerspring and slat frame beds. Installation is quick and easy and the LegLifter can be specified for either the left or right hand side of the bed.



## Specifications

LegLifter unit dimensions	Width 1050mm Height 660mm
Width of leg-rest	485mm wide
Height adjustment range	500mm to 700mm
Total weight including packaging for transit	18kg
Weight of LegLifter unit	9kg
Installation time for standard domestic bed	From 5 to 10 minutes
Maximum leg-load weight	Tested to 67kg: legs must fit within the width of the leg-rest
Distance from bed-frame to platform of folded leg-rest	60mm

Order codes	
Left side of bed	MK3L - LH
Right side of bed	MK3L - RH

\* Ask for the detailed brochure on this product



# The Multifit Guarantee

---

## Covers all our chair and bed raisers

Safe and secure.

Long lasting.

Robust, re-usable, economical.

Telescopic adjust-ability for perfect fit.

Unique Multifit Tri-knob will not loosen.

Delivery within 5 working days from order or we will treat your entire department to a morning tea.

The equipment will fit or we will customize one that will.

The equipment will retain its original function for at least two years or we will replace it at no extra cost.

Quality of life is precious. Our mission is to provide you with quality long lasting equipment that enables you to live your life with independence.