

Bed Blocks




Fig 1

- (1) Perimeter lip, (2) Recess 75 x 75 wide x 40mm deep
 (3) Finger hole, (4) Pyramid shape

Raising the height of your bed may enable you to get in and out of the bed more easily.

- These light weight bed blocks have finger holes to make the blocks easier to carry
- Recess has a non slip rubber pad
- Durable construction means they can be reused over and over again while maintaining excellent condition
- Pyramid shape means that they are self-bracing
- Sealed surfaces mean cleaning is made easier


- Re Issuable
- Light weight
- Cost effective
- Stable pyramid shape

SPECIFICATIONS

Max user	380Kg	
Perimeter lip	40mm	
Pad size	75 x 75 mm	
Quantity	Simply order the number of Blocks to match the bed	
Bed Size	Will fit any bed where the bed legs are < 75 x 75 mm	
Bed Blocks		
Height	Supplier Code	
150mm	MF 0/4	
100mm	MF 0/3	
75mm	MF 0/2	
50mm	MF 0/1	

Guarantee

- You will receive any list item within five working days or we will treat you morning tea.
- Our equipment will fit or we will customise one for you that will.
- Our equipment will retain its original function for at least two years or we will send you a replacement at no extra cost.



Pg1



Multifit

Get in & out of your bed with less effort.



Fig 2

Bed Blocks will fit this type of bed




Fig 3

Bed Blocks fitted



Fig 4

If you require increased stability choose the Multifit Bed Raiser where the Blocks are linked with a telescopic section.



LARGER FEET?

Choose the Widened Bed Raiser

Fig 5



SMALLER FEET ?

Choose the Divan Bed Raiser

Fig 6




Made with pride
 In New Zealand
 for over 25 years



Multifit Hospital Supplies Ltd © 2020,
 Ph: 0800 102 803, E: reception@multifit.co.nz
 Website: www.multifit.co.nz 250720

Bed Blocks

Fitting Instructions



MultiFit

Get in & out of your bed with less effort.

The following method ensures Bed Blocks are fitted correctly, safely and with the minimum of effort.




Using the Multifit Bed Jack (1) lift the end of the bed up. Fit the Bed Blocks (2) centrally under the feet of the bed. Lower the Bed Jack.



Place a weight (3) to prevent the bed tipping. With the Multifit Bed Jack (4) Lift the head end of the bed. Fit the Bed Blocks (5) centrally under the feet of the bed. Lower the Bed Jack.



Confirm the raised bed is safe and stable.



Check:
Is getting in and out of the bed easier?
Was the bed raised to the correct height?
Does the user look confident?

Maintenance - check monthly the stability of the raised Bed. Repeat steps A,B,C,D if required.

Cleaning - with damp cloth and non-abrasive detergent, clean and dry thoroughly.

Demonstrations - see our stand at the next Show Your Ability Expo.

Parts replacement - refer to www.multifit.co.nz for a complete parts list.

Our Mission - is to provide you with the safest, most cost-effective raiser to meet your needs.



Scan for more info



Made with pride
In New Zealand
for over 25 years



Multifit Hospital Supplies Ltd © 2020
Ph: 0800 102 803, E: reception@multifit.co.nz
Website: www.multifit.co.nz 250720