

5 Leg Chair Raiser MF 1/Q




Fig 1

The 5 legs of the raiser adjust to fit the legs on the chair. Fitting Instructions.

Simply extend the raiser in width to suit the frame on the base of your chair (1- fig 2)

Adjust the Raiser so that it fits tightly up to the frame of your chair (2 - fig 3)

Tighten the 5 x clamp bolts on the Raiser (6 - fig 3)

Loop the cord over the leg of the chair (3 - fig 3)

Thread the cord through the cleat (4 - fig 3)

Tighten and fix the cord in the cleat (5 - fig 3)

Finally Test for safety and stability before use.



Raiser will fit this type of chair

Specifications

Max user weight	380Kg (+ the furniture)			
Pad Size	90 mm square			
lip around the perimeter	20mm			
5 Leg Raiser Chair Platform				
Height	Size	Supplier Code	Order Code	JDE#
150mm	650 – 820 mm	MF 1/Q-150		
100mm		MF 1/Q-100		
75mm		MF 1/Q-75		
50mm		MF 1/Q-50		
Optional: non slip rubber pads stop raiser slipping on polished timber floors				
Safety First - Raiser must match the chair type. See www.multifit.co.nz				
Tested - to ensure the stability of the raised furniture. This safeguards the user.				




Fig 3

Expand the raiser out and fit it to the base of the chair.

Assessment notes: Was the chair raised to the correct height? Does the user look more confident? Chairs are most stable when getting in and out. Reaching over the side of the chair is not recommended. If over time the situation changes and the user needs additional assistance, consider Chair Transfer Rails or an orthopaedic chair.

Maintenance - check tightness of the raiser to the chair, tightness of wing nuts, knobs and fixings.

Clean - with damp cloth. Use non-abrasive detergent and dry thoroughly.

Demonstrations - see our stand at the next Show Your Ability Expo.

Parts replacement - refer to www.multifit.co.nz for a complete parts list.

Our Mission - is to provide you with the safest, most cost-effective raiser to meet your needs.



Made with pride
In New Zealand
for over 25 years



MultiFit Hospital Supplies Ltd © 2020
Ph: 0800 102 803, E: reception@multifit.co.nz
Website: www.multifit.co.nz 180720