



AquayJoy Bathlift—Note that the back rest and therefore the child needs to be upright when the Bathlift is raised and that the wings on the Bathlift are fixed not hinged. That means that the child must be placed onto the seat then moved back between the wings. Then the you fix the pommel. And the straps.



Mangar Surfer Bather— note that the back rest is reclined when the Bathlift is raised and lowered. The back rest is curved, the horse shoe shaped head rest and the lap strap are adjustable in height. All of these things help to hold the child centrally on the Bathlift and make it easier on the care giver to place the child on and off the Bathlift. There is a limitation on the size of child that the Bathlift suits though and this is determined by the length of the childs legs and where the taps are located in the bath. EG: As the Bathlift raises up care must be taken to ensure the child's feet do not contact the taps. Max user height is 1350mm and 50kg