

Safe-er-Grip Technology

Install Instructions and Trouble Shooting:

Install:

1. Safe-er-Grip products may only be applied to smooth, non-porous surfaces. No texture or painted surfaces of any kind can be present.
2. Prior to installing all surfaces must be cleaned thoroughly.
3. If applying to tile, tile squares must be 4" x 4" or greater.
4. Flip Up release latches.
5. Position Safe-er-Grip where desired on tub or shower wall.
6. Suction cups must not overlap grout lines.
7. Apply pressure to the handle and flip latches down to secure onto wall.

Remove:

1. Release both latches.
2. Slide fingertip under rubber seal at each end if necessary, to release suction.

Maintenance:

1. Must re-attach Safe-er-Grip product prior to use each and every time.
2. Always test suction on Safe-er-Grip product prior to use.
3. Make sure surfaces underneath suction cups are cleaned regularly.

Trouble Shooting:

1. Product will not attach to tub or shower wall:
 - a. Make sure tub or shower wall is a smooth, non-porous surface. This product will not work on textured or painted surfaces.
 - b. Make sure suction cups on Safe-er-Grip product are not covering any grout lines. No air can be allowed to pass under suction cup or it will not hold.
 - c. Make sure surface is thoroughly cleaned prior to installation.