

● Accessories

- Arms with or without an integral operating switch offer additional support. They tilt up for easy side transfer or for others to use the toilet.
- Extended feet allow the lift to be placed over higher toilets (additional 64mm/2.5" height)*.



* Must be fitted by Moutway when ordering the lift.

Technical Specifications

Maximum user weight	160kg/25 stone/350lbs
Total product weight (incl. battery but no arms)	19.04kg/42lbs
Overall width without arms	612mm/24.09"
Overall width with arms	710mm/27.95"
Width between arms	530mm/20.8"
Seat height (down)	430mm - 490mm /16.9"-19.3"
Seat height (up)	560mm - 620mm/22.05"-24.41"
Toilet bowl clearance height	370mm - 430mm/ 14.56"-16.9"

Contact your local retailer/supplier or visit www.moutway.com for full product specifications.

● Order codes

Item Order	Code
Toilet lift with no arms (hand control operation)	MS0HSUK
Toilet lift with 2 arms (left operating switch)	MS2LSUK
Toilet lift with 2 arms (right operating switch)	MS2RSUK
Left / right hand arm (no operating switch)	TL626/S (left) TL627/S (right)
Left / right hand arm with operating switch	TL625/S (left) TL628/S (right)
Extended feet	TL634/S

● Guarantee



The Moutway Solo Toilet Lift is guaranteed for 2 years and the batteries for 1 year.

For more information or to order:

call: **0800 014 1852**
e-mail: mail@moutway.co.uk
or visit: www.moutway.com



With 30 years' experience, Moutway prides itself on 'bringing independence to life' with innovative products such as the Moutway Solo Toilet Lift.

Visit www.moutway.com or ask your local retailer/supplier for information on our other product solutions for independent living.



moutway
solo[®]
toilet lift

Replicates natural standing motion

Modern, stylish design

No complicated installation

Maintain your independence and dignity in style!



bathing and toileting



Part number: ML1042 / Issue: 1; 03/11

The Mountway Solo Toilet Lift is ideal if you struggle onto or off the toilet, perhaps due to muscle weakness or stiff joints.

At the simple touch of a button, the lift will lower you gently onto the toilet, then bring you effortlessly back up into a natural standing position.

If you have carers present, they won't need to lift, as the powered seat does all the work!



● Modern, stylish design

With its smooth lines and easy-to-clean surfaces, the Mountway Solo Toilet Lift blends discreetly into any bathroom setting. In its lowest position, the lift becomes a standard toilet seat easily used by others in your household.

● Replicates natural standing motion



A sturdy seat provides excellent support to safely lower you onto the toilet. It then brings you smoothly upwards and forwards into a standing position. For your optimum comfort, the seat can be stopped in any position*.

** Standard operation by a lightweight, easy-to-use hand control.*

● No complicated installation

Simply remove your existing toilet seat and place the lift over the toilet bowl for immediate use. Adjustable feet allow positioning on uneven flooring and over toilets of various heights.



● Safe and simple to use

The Mountway Solo Toilet Lift is completely battery-powered and is supplied with two rechargeable batteries for uninterrupted use. A built-in battery monitoring system will also tell you when recharging is needed*.

** Battery charger also supplied*

