

Bed Rail Std

MF 3/1 STD

Used to assist in sitting up or turning in bed as well as transferring in and out of bed. This Rail can also be used to pull on to assist in standing if you have fallen beside your bed.

- 16mm base board is specially shaped to give you maximum stability.
- The rail is welded to the base board for superior strength.
- All the edges are rounded to protect you from injury
- Welded joints enable the bed rail to fit to standard beds as well as slat beds.
- Available with a cord and cleat or strap and buckle to prevent the Rail pulling out from the bed. This feature gives the bed rail a very solid feel and the security of knowing you are safe.
- Powder coated rail and clear coated base board give the rail a long maintenance free life

Maximum user weight of 150kg makes it one of the strongest rails on the market and ensures your safety when using it.



Technical specifications :

Maximum user weight	150kg
Weight	7 kg
Rail diameter	25 mm steel tube
Height of rail	50 cm or 60 cm options
Width	46 cm
Length	79 cm
Guarantee :	2 years



Can be used to pull on to assist in standing if you have fallen beside your bed



All the edges are rounded to protect you from injury

Max user weight of Bed Rail 150 kg is hard to beat.



Whenever you fit a Bed Rail its a good idea to have it strapped to the bed to stop it from pulling out. Normally the Bed Rail is used to assist in sitting up or turning in bed as well as transferring in and out of bed. But occasionally people fall beside the bed and use the Bed Rail to pull on to assist in standing. Used in this way the Bed Rail cannot pull out from the bed, but a huge amount of load is transferred onto the Bed Rail and its important to know that it will be strong enough.

We recently tested our Bed Rail by imposing a load equivalent to a person of 150kg by pulling on the rail at an angle of 45 degrees. Our Bed Rail came through with flying colours proving yet again that the Multifit Bed Rail is one of the strongest on the market today.