## **Bed Lever MF 3/1AS-CC**










Placed between the mattress and the base of your bed, the Bed Lever provides a useful handhold to assist you on and off your bed safely. Correct fitting is essential to the safe use of this product. Take care to read the fitting instructions and ensure that it is refitted the same way every time the bed is made.

Max user weight	150kg (25 Stone) when fitted to a double bed or larger	
Loop size	550 H x 225W (81mm between the vertical tubes)	
Base size	465 W x 835 L x16 D (all measurements in mm)	
Suits bed types	Divan & Inner-sprung beds (all sizes). Does not suit profile beds	
Cleaning	Clean with a soft cloth moistened with soapy water.	
		MF 3/1AS-CC
		Cord & Cleat
	be released by pulling on the cord	Cord & Cleat
where it extends	under the handle.	Cord & Cleat
where it extends This means the Bo	under the handle.	Cord & Cleat

## Keep the Loop handle away from the pillow



Phone: 0800 102 803, Email: reception@multifit.co.nz Website: www.multifit.co.nz