

Bed Jack

Code : MF BJ / BAT (Mk3)



Reduce Manual Handling risks when fitting bed raisers.

Simply place the foot of the jack under the edge of the bed and push the button to raise or lower your bed.

Latest improvements are : Wheels for better portability, lengthened lift plate so it reaches further under the bed, new battery housing and switch making them both more robust.

- Powerful lifting action lifts any bed with the touch of a button.
- Battery will last all day on a full charge.
- All moving parts are guarded to improve your safety.
- Handle on top is used to place the jack.
- Battery powered means you don't have to worry about leads
- Wheels on the bottom make it easily portable.
- Narrow profile means it can be placed in tight spaces .

For a detailed New Zealand case study on the benefits of using a Bed Jack to lift a bed when fitting bed raisers, go to our website www.multifit.co.nz



Simply push the button to lift the bed



Battery charger simply connects to the top of the battery and charges over night.

Specifications :

Max load it will lift	250 kg
Lift stroke	20 - 600 mm
Frame Colour	White
Width	210 mm
Weight	7 kg
Height	870 mm
Construction	Light weight steel tube and aluminum extrusion
Battery	
Recharge	Trickle charge over night
Lifts per charge	50 lifts -300mm stroke 35 Kg load

Multifit
Hospital Supplies Ltd

Multifit Hospital Supplies Ltd
10 Pohutu Street,
Whakatane
New Zealand.
Phone 0800 102 803
Fax 0800 102 800
E mail : info@multifit.co.nz
Web site : www.multifit.co.nz

©2005 Multifit Hospital Supplies Ltd